



QuikRead go[®] Instrument software 7.6.2 release note

A new QuikRead go Instrument software version 7.6.2 will be released in January 2019. With the new software version it will be possible to measure the new QuikRead go[®] easy CRP test.

The new software version has been developed in compliance with IEC 62304:2006+AMD1:2015 Medical device software -- Software life cycle processes and General Principles of Software. It has been validated according to Final Guidance for Industry and FDA Staff Jan 11 2002. Testing included measurement of 33 samples with both software versions, new 7.6.2 and previous 7.5.1. All results passed.

An instrument with any of the previous software versions can be updated to this new software version. Version update does not delete old results from the memory of the instrument.

Features

New analytes

- QuikRead go[®] easy CRP

Improvements

- More quiet mixing sound with QuikRead go[®] wrCRP

Availability of new software

- You can download the new software version 7.6.2 from <http://www.softwareupdate.quikread.com/>
- The new software version can be also ordered on a USB memory stick, cat.no 137540 QuikRead go Software USB stck

